



Cleaning and Disinfection of Respironics Profile™ Lite Mask

CLEANING (Perform prior to proceeding to any of the disinfection or sterilization methods listed below.)

- ✓ Remove the headgear and spacer, if appropriate.
- ✓ Soak the parts in a commercially available enzymatic cleaner (e.g., Klenzyme™).
- ✓ Clean the Profile Lite Mask with a soft bristle brush in a solution of cool tap water and a commercially available anionic detergent (e.g., Manu Klenz™). Do not use cleaning products that contain conditioners or moisturizers because they will leave a residue.
- ✓ Rinse thoroughly under cool running tap water. Air dry.
- ✓ If adhering substances are not adequately removed, replace the mask.
- ✓ See options below for disinfecting or sterilizing the Profile Lite Mask.

PASTEURIZATION

Temperature: 70°C
Exposure time: 30 minutes
Result: High level disinfection
No. cycles tested: 30

- If the swivel does not move freely after disinfection, replace the split washer.*
- Inspect the mask after processing. If any components are damaged, replace the Profile Lite Mask.

CIDEX® (Do Not Use CIDEX O.P.A.)

Exposure time: 45 minutes
Rinse: Per manufacturer's instructions
Result: High level disinfection
No. cycles tested: 30

- If the swivel does not move freely after disinfection, replace the split washer.*
- Inspect the mask after processing. If any components are damaged, replace the Profile Lite Mask.
- Slight discoloration of the mask cushion after processing is normal.

STERRAD® System 100

Exposure time: Machine controlled
Result: Sterilization
No. cycles tested: 30

- If the swivel does not move freely after sterilization, replace the split washer.*
- Inspect the mask after processing. If any components are damaged, replace the Profile Lite Mask.
- A slight odor may be present after processing. This is normal.

*Part #1002720 (package of 5 washers)

Profile is a trademark of Respironics, Inc. CIDEX and STERRAD System are registered trademarks of Advanced Sterilization Products.

- Inspect the mask after processing. If any components are damaged, replace the Profile Lite Mask.
- Klenzyme and Manu Klenz are trademarks of Steris Corporation.*